

Vario 18 c/i



RELIABLE CONSTANT OR INTERMITTENT
SUCTION



Vario 18 c/i Suction Pump

Vario 18 c/i

Reliable constant or intermittent suction

The Vario 18 c/i is a high-quality medical suction pump for constant or intermittent use. It is powered by the well-proven QuatroFlex system that allows for easy servicing, plus it combines easy handling and cleaning with safety features to ensure optimal operation. The Vario 18 c/i guarantees maximum suction performance for many constant or intermittent suctioning needs, such as gastric suctioning. Intermittent suction can provide a more efficient way of drainage at times and is ideal for delicate procedures, whereas the constant suction mode can be used for airway suctioning or smaller interventions. The versatility of the device makes the Vario 18 c/i an all-rounder for small surgeries, airway suctioning and intermittent use for hospitals and clinics.

Vario 18 c/i can be used with either a disposable or a reusable collection system, according to each customer's needs and demands.

Manufactured with Swiss precision and quality, Vario 18 c/i offers a range of useful features:



- a** Overflow protection/bacteria filter
- b** Vacuum gauge in kPa and mmHg
- c** LED warning light
- d** Carrying handle
- e** Membrane vacuum regulator
- f** Tubing holder
- g** Plug for 12 VDC IN car connection cable
- h** Constant/intermittent/OFF switch

Distinctly Swiss:

The innovative 4-fold piston/cylinder system, QuatroFlex, builds up the vacuum within seconds. As well as being quiet in operation and unique in design and technology, it is also easy to service and maintain which is highly appreciated by medical technicians. The membrane vacuum regulator allows for precise vacuum settings. It needs to be pressed down to change the vacuum, thus making it impossible to change the settings inadvertently. The Vario 18 c/i can therefore be used at the patient's bedside.



Precise vacuum setting with press-down knob

Versatility for a wide range of applications

The Vario 18 c/i is lightweight and portable and can be used in many departments in hospitals and clinics. Variable suction levels can be set with the vacuum regulator so the Vario 18 c/i can be used for surgical interventions or for suctioning procedures that require intermittent settings. With the battery version, the Vario 18 c/i is mains-independent for 30 minutes. It can also be equipped with either a disposable or a reusable Medela collection system.



Can be equipped with either a disposable or reusable collection system

Ease of use

Vario 18 c/i offers intuitive design and easy handling. All accessories are designed to complement it perfectly and to prevent false connections. The Vario 18 c/i is even more flexible with optional accessories including a trolley, carrying bag and various attachments. It features an inbuilt battery charger that allows for mobile use. Acoustic and visible alarms notify the user when the battery needs recharging.



Optimal hygiene with overflow protection/bacteria/odour filter

Fast Facts

- Our Swiss development and manufacturing process yields mechanical precision to guarantee excellent performance.
- QuatroFlex technology: silent function: can be used in environments in which quiet operation is necessary.
- Precise vacuum settings in the medium vacuum range and fast vacuum build-up: ready when needed.
- One device for multiple needs.
- Mains-independent for 30 minutes.
- Acoustic and optical battery management.
- Wide range of accessories available for individual application.
- Vario 18 c/i incorporates reliability and ease of use like the Medela suction pumps for operating theatre.
- 2 years warranty on device.

Delivering excellence for over 50 years with comprehensive, high quality Swiss medical vacuum systems.

Accessories

Filters

to protect the pump from overflow, the environment from bacteria and to neutralise odour.



Trolley

with rails for mobile transport of the pump and accessories.



Reusable Collection System

autoclavable jars and lids collect secretions easily and economically. Sizes available: 1 and 2 litres.



Holders and plates

can be used to attach Vario 18 c/i to standard rails.



Disposable Collection System

with disposable suction liners and reusable PC jars that offer a hygienic, simple handling and efficient fluid management. Sizes available: 1.5 and 2.5 litres.



Carrying bag

for transporting the pump.



Use Medela accessories for even greater flexibility and ease of use.

Plug for 12 VDC IN car connection cable

for operating the pump in the car.



Technical data

	Medium vacuum -55 kPa - 413 mmHg (Tolerance: +/- 10 %)
	low flow 18 l/min. (+/- 10 %)

	AC 3,5 kg AC/DC 4,2 kg
	H x W x D 380 x 170 x 285 mm

	IP21
	IEC 60601-1:2005, Edition 3.0 fulfilled 2 years warranty
	

Medical Vacuum Technology for Healthcare Professionals

Please contact us or your local Medela representative for details.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Local contact:

Medela Australia Pty Ltd
Medical Technology
3 Arco Lane
Heatherton, Vic 3202
Australia
Phone +61 (0) 3 9552 8600
Fax +61 (0) 3 9552 8699
contact@medela.com.au
www.medela.com.au

Medela Canada Inc.
4160 Sladerview Crescent
Unit #8
Mississauga, Ontario
Canada L5L 0A1
Toll Free 800 435 8316
Toll Free Fax 800 995 7867
info@medela.ca
www.medela.ca

Medela UK Ltd.
Huntsman Drive
Northbank Industrial Park
Irlam, Manchester M44 5EG
United Kingdom
Phone +44 870 950 5994
Fax +44 870 389 2233
info@medela.co.uk
www.medela.co.uk

Medela Inc. USA
1101 Corporate Drive
McHenry
IL 60050
USA
Phone +1 877 735 1626
Fax +1 815 363 2487
suction@medela.com
www.medelasuction.com

Medela India Office
Plot number 52,
Udyog Vihar Ph-VI, Sector-37,
Near Hero Honda crossing,
Gurgaon-122001
Tel. 0124 - 4121600
Toll Free: 1800-102-3645